ADAPTACIÓN LENTES MULTIFOCALES

MULTIFOCAL CONTACT LENS FITTING

LENTE MULTIFOCAL SERVILENS

SERVILENS MULTIFOCAL CONTACT LENS

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Simultaneous Vision Design
Art or Technology

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Simultaneous Vision
Elderly or Young

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
Servilens Multifocal Lenses Success

Reproducibility of Multifocal Lenses.

Servilens Multifocal Software was running automatically during one week to calculate empirically more than 500.000 different Multifocal Contact Lenses and to register a complete range of Multifocal lenses for the eye care professional.

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Three Add Ranges to fit more Presbyopes.

A complete range of Soft toric multifocal lenses. (4 days)

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MULTIFOCAL CONTACT LENS FITTING

Suggestions

-Indicate to the patient the differences between the selected multifocal contact lens and spectacles.

-First comfort and then vision. The contact lenses should be worn for at least 4 continuous hours straight and the patient should be asked to identify any comfort problem which might be occurring related to the contact lens wear (biomicroscopic examination).

-Make use of proper illumination when carrying out visual tasks.

-Patient with astigmatism (greater than 0.75 diopters) should wear multifocal toric lens because acuity is critical in multifocal contact lens fitting. NIPO test is very important for multifocal toric lens fitting.

-In some cases, having supplemental spectacles to wear over the multifocal contact lenses for specific visual tasks may improve the success of the fitting (i.e. when driving at night).

**STEPS TO FITTING:
First comfort and then vision**

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Multifocal contact lens fitting.

1.- Anamnesis

2.-Determine which is the dominant eye prior to a new refraction:

Method 1= Fogging with a +2.00 D lens placed alternately in front of the right and left eye. Dominant eye will be the eye with more distance visual annoyance.

Method 2: Have the patient point to an object at the far end of the room and cover one eye. If the patient is still pointing directly at the object, the eye being used is the dominant eye.

3.- New refraction for fitting CL.

Far vision refraction: considering dominant eye for far vision.

Far vision refraction: considering non dominant eye for near vision.

Near vision refraction: (working distance) when possible, choose the minimum Add.

Case Hyperopia:

We recommend a +3.00 fogging test (and reduce powers) to prescribe maximum plus power for Far Vision.

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Case Myopia:

We recommend to prescribe minimum power for far vision.

- Case astigmatism: confirm the best axis (5° or 10°)

4.-Pupil diameter in normal light, Corneal evaluation, tear film analysis.

5.-Corneal measurements must be taken :

6.-Order the contact lens, correcting for the back vertex (far vision). Order CL online at: www.lens55.com

Design and Add power determination:

Initial dispensing:

Servilens recommends to fit the same multifocal design for both eyes to improve binocular vision.

Design 1 (Dominant Center Distance)

Patient Visual Demands
Far vision more important than Near vision. e.g. driver.

Design 2 (Non Dominant Center Near)

Patient Visual Demands
Near vision more important than far vision. e.g. secretary.

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Soft Toric Multifocal Over-refraction

Over refraction of: +0.50 or +0.75 on spherical multifocal lens (far or near)

You can add to the final CL power.

Over refraction of: +0.50 or +0.75 on toric multifocal lens (far or near)

You can not add directly to the final CL power because:

- It could be due to spherical optical surface (if you add it you are right).

or

-It could be due to misalignment of the cylinder axis optical surface (if you add it you could be wrong). The problem is the axis. Apply LARS

Recommended solution

First NIPO test (see lens55.com) and apply spherical over-refraction (+0.50 or +0.75) with two cylinders fixed in the trial frame.

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MULTIFOCAL CONTACT LENS FITTING

Special Fitting Other Suggestions

Anisometric patient (~3.00)

More myopic eye for near vision.

Emmetropic patient (Add 1.75)

Possible latent hyperopia.

If not: fit only one contact lens on the non dominant eye (+1.75).

Patient (RE -2.50 D / LE -1.50 D Add 1.75)

Fit one only contact lens on the right eye (-2.50).

When critical vision is needed at a certain distance, fit a monofocal contact lens to the required distance.

IFU Instructions For Use- Contact
Lens Care

Servilens recommends to deliver and sign the IFU. Once the right contact lens has been chosen, instructions for lens insertion and removal, hygiene, and wearing time will be provided for the user.

ADAPTACIÓN MULTIFOCALES

MULTIFOCAL FITTING