ADAPTACIÓN RGP

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B- Initial Lens Power Selection.

1) Convert Rx to minus cylinder if necessary.

2) Correct for vertex distance if either meridian is greater than ± 4.00 using vertex distance chart.

3) Power will be equal to spherical component of the spectacle correction (corrected for vertex distance expressed in minus cylinder format) "On-K" fit. Corneal astigmatism compensation 3.50 D.

4) RULE=SAMFAP (Steeper,Add Minus Flatter Add Plus) correction must be made if the lens is steeper or flatter than K (or than the trial lens used). Change the power by the dioptric equivalent of the change in base curve.

Characteristics of a Well-Fit Spherical/Aspheric Lens

- The lens should center well over the pupillary zone on the cornea.
- The lens should move freely with the blink.
- The fluorescein pattern should show good tear exchange.

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
RGP FITTING

- After wearing 8 hours, instill fluorescein into the eyes and exam a health cornea

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

RGP FITTING

KERATONOCUS-POST LASIK FITTING

The central keratometer measurement does not help in some cases, especially in post-refractive surgery and keratoconus (and some hyperopes). A corneal topography map provides more information. Eccentricity is very important. The choice of the first diagnostic RGP lens in these cases should be the average of the superior hemisphere corneal values. We would recommend to fit a loose lens and later on, if necessary, to try a tighter Base Curve.

Remarks:

Fit flatter lenses whenever comfortable, good vision and the corneal epithelium is perfect after 8 hours of continuous use.

The best choice is fitting the RGP lens "ST K 40 (Keratocoonus)" together with "Lens 55 Silicone 1 Day " as PIGGY-BACK.

Other option: Lens 180 GP+S. Hybrid lenses.

RGP
Fitting